

# TPOL Hardware Fall 2019 Checkup



# TPOL supplies

Nick found all of the TPOL supplies and consolidated the items in one location

## Spares

- Three Swan preamp boxes
- Three functional filter boards (2 from ASU 1 from JLab)
- Two blank filter boards
- 1 blank board for distribution box

# Location of TPOL supplies



**Tool  
box**



**Cabinet**




**TPOL supplies at  
bottom of Cabinet**

# Signal check



Checked each channel for signal

- Typical result shown below
- All channels included at end of this presentation





# Baseline and noise check

- Baseline fluctuation is typically a few mV and noise envelope about +/- 5 mV
- Typical example shown below (all at end of this presentation)
- Worst outliers shown on next slide






# Baseline and noise check

- Channel 10 noise envelope is about twice the size as compared to a typical channel
- Baseline for channel 32 is a very clear periodic fluctuation of about +/- 5 mV





# Noisy Channel 10, 5s envelope, 2.8mV trig



Zoom in on envelope shows that time structure is not very sharp (voltage versus time has low slope)

# Noisy channel 10, 30 mV trigger



- Signal has sharp time structure



# Trigger


Conversions:

- 5.33 keV per mV
- 2.6 keV per FADC channel

A threshold of 240 keV implies 45 mV from bottom of signal to trigger


A trigger set to 30 mV would require a signal bottom within 15 mV of nominal

# Noisy channel 10, 30 mV trigger




- Want distance between baseline and bottom  $< 15$  mV

# Large baseline channel 32, 30 mV trigger



- Want distance between baseline and bottom  $< 15$  mV

# Typical channel 3, 30 mV trigger



- Need distance between baseline and bottom < 15 mV





# Remarks

- Did not find any need to switch out the preamps or filter boards 😊
- All looks well for TPOL to run using a 30 mV trigger





All scope pictures shown in remaining  
slides



# Channels 1-4







# Channels 5-8













# Channels 9-12







# Channels 13-16






# Channels 17-20







# Channels 21-24







# Channels 25-28







# Channels 29-32



# Channels 1-2







# Channels 3-4













# Channels 5-6







# Channels 7-8







# Channels 9-10







# Channels 11-12







# Channels 13-14







# Channels 15-16



# Channels 17-18







# Channels 19-20













# Channels 21-22







# Channels 23-24







# Channels 25-26







# Channels 27-28



# Channels 29-30



# Channels 31-32



# Title

